




The Gentle Choice for Accurate Therapy

The Bubble CPAP System is a convenient method to deliver Continuous Positive Airway Pressure (CPAP) to neonates and infants requiring respiratory support.

Used alongside the FL-9000 Respiratory Humidifier Heater, the Bubble CPAP System provides essential humidification during the whole therapy.

The extensive size range of Flexicare's Infant Nasal CPAP prongs provide optimal care for the most vulnerable patients.


Bubble CPAP System



038-33-180	Bubble CPAP System - KIT 1 Neonatal 15mm Inspiratory Limb Heated Wire 1.6m, 0.3m Incubator Extension, 10mm Expiratory Limb 1.2m, Bubble CPAP Generator, Pressure Relief Valve 40 cmH ₂ O, Universal 22mm Connector, 15F/8M Connector, 15F/10M Connector And Auto-Fill Chamber.	10
038-10-300	Bubble CPAP Generator	10

Infant Nasal CPAP Interface

Infant Nasal CPAP includes colour-coded Nasal Cannula, Soft Knitted Hat with securing straps and Sizing Guide.



PART NUMBER	DESCRIPTION	SIZE	APPROX WEIGHT RANGE	COLOUR	QUANTITY
038-13-100	Infant Nasal CPAP	Size 00	Less than 500g	Blue	10
038-13-110	Infant Nasal CPAP	Size 0	500 – 700g	Black	10
038-13-111	Infant Nasal CPAP	Size 1	700 – 1250g	White	10
038-13-112	Infant Nasal CPAP	Size 2	1250 – 2000g	Green	10
038-13-113	Infant Nasal CPAP	Size 3	2000 – 3000g	Yellow	10
038-13-114	Infant Nasal CPAP	Size 4	Over 3000g	Red	10
038-13-115	Infant Nasal CPAP	Size 5	1 to 2 years old	Orange	10