



Clic™ Closed Suction Catheters


Protecting patients and caregivers

Effectively clears patient airways of accumulated mucosal secretions; improving respiratory function, integrated Bacterial/Viral filters protects clinicians from air dispersing at the Clic thumb valve.




Reduce Contamination Risk

Traditional Closed Suction Catheters leak contaminated air during the insertion of catheter into the patient airway, allowing viruses and bacteria to escape into the clinical environment.



Added layer of Protection

Clic features 99.99% efficiency Bacterial/Viral filters, creating a secure barrier by filtering the escaping contaminated air, protecting both clinicians and patients.




PRN Up To 7 Days

Extended use of Clic, reduces cost and creates less waste.



Optional Ports

Optional MDI port for medication delivery and EtCO₂ port for carbon dioxide monitoring. In both T-piece and Double Swivel Elbow configurations.



Comprehensive Range Of Closed Suction Catheters

All sizes and configurations to meet all patient population needs.

SCAN HERE



3D VIEW

Clic Closed Suction Catheters maintain a continuous closed system, protects ventilation and pressure levels, and reduces risk of cross-contamination and Ventilator Associated Pneumonia (VAP).

Clic Range of Closed Suction Catheters:

WITH T-PIECE	WITH DOUBLE SWIVEL ELBOW	WITH DOUBLE SWIVEL ELBOW AND MEDICATION PORT	WITH DOUBLE SWIVEL ELBOW AND ET/CO ₂ PORT	DESCRIPTION	QTY
050-63-110	050-64-110	050-65-410	050-65-310	Adult ET Closed Suction Catheter Size 10Fr 540mm (21.3")	20
050-63-112	050-64-112	050-65-412	050-65-312	Adult ET Closed Suction Catheter Size 12Fr 540mm (21.3")	20
050-63-114	050-64-114	050-65-414	050-65-314	Adult ET Closed Suction Catheter Size 14Fr 540mm (21.3")	20
050-63-116	050-64-116	050-65-416	050-65-316	Adult ET Closed Suction Catheter Size 16Fr 540mm (21.3")	20
Y-SETS				DESCRIPTION	QTY
050-62-105				Neonatal Closed Suction Catheter Size 5Fr 305mm (12") with Y-Sets 2.0/2.5mm	20
050-62-106				Infant Closed Suction Catheter Size 6Fr 305mm (12") with Y-Sets 2.5/3.0/3.5mm	20
050-62-107				Infant Closed Suction Catheter Size 7Fr 305mm (12") with Y-Sets 3.0/3.5/4.0mm	20
050-62-108				Infant Closed Suction Catheter Size 8Fr 305mm (12") with Y-Sets 3.5/4.0/4.5mm	20
050-62-110				Paediatric Closed Suction Catheter Size 10Fr 405mm (16") with Y-Sets 4.5/5.0/5.5mm	20
WITH T-PIECE	WITH DOUBLE SWIVEL ELBOW			DESCRIPTION	QTY
050-63-210	050-64-210			Tracheostomy Closed Suction Catheter Size 10Fr 305mm (12")	20
050-63-212	050-64-212			Adult Tracheostomy Closed Suction Catheter Size 12Fr 305mm (12")	20
050-63-214	050-64-214			Adult Tracheostomy Closed Suction Catheter Size 14Fr 305mm (12")	20
050-63-216	050-64-216			Adult Tracheostomy Closed Suction Catheter Size 16Fr 305mm (12")	20